



SAVE THIS MANUAL: KEEP THIS MANUAL FOR SAFETY WARNINGS;  
PRECAUTION OPERATING, INSEPTION, MAINTENANCE AND CLEANING PROCEDURES.

# IMPORTANT SAFETY INFORMATION

## WARNING

### GENERAL SAFETY WARNINGS

Read all safety warnings and instructions. Failure to follow the warnings and instructions may result in electric shock, fire and/or serious injury. Save all warnings and instructions for future reference. The warnings, precautions, and instructions in this manual can not cover all possible conditions and situations that may occur. It must be understood by the operator that common sense and caution are factors which cannot be built into this product and must be supplied by the operator. Read carefully and understand all assemble and operation instructions before use of this product. Failure to follow safety rules and other basic precautions may result in serious personal injury.

- For adults only.
- Use the massage gun on dry and clean surface of the body. Press and move gently on the skin.
- Massage time on the same muscle parts should not exceed one minute at a time.
- Do not use the massage gun on your head or bony parts.
- Choose a suitable massage head for the specific need. Frequent high-speed pressure on one place can cause abrasions. Stop use immediately if you experience any pain during massage.
- Keep fingers and hair away from the massager to prevent injury.
- Keep massager away from heat sources.
- Keep the massage gun away from water.
- Inspect the massage gun for any damage prior to each use. If damaged, do not use the massager.
- Modifications of this massage gun is strictly prohibited.
- Do not use high-speed gear / equipment without medical approve if: 1. Pregnancy, Diabetes and complications, no neuropathy or optic lesions, using a pacemaker, have had surgery, have epilepsy or migraines, disc herniation, relaxation or spondylitis. 2. Recent joint replacement or IUD, or physical problems. Weak adults and children need to use this machine carefully, accompanied by a healthy adult.
- Contact a physician if you are unsure about using the massage gun.

# OWNER'S MANUAL AND SAFETY INSTRUCTIONS

## MESSAGE GUN PRO



Motor Type: Brush Motor  
Charging voltage: 5V Adaptive Power  
Regulation Voltage: 7.4V  
Rated power: 25W; 6 gears Motor speed: 1800-3200  
Pool: 1500\*2/1800\*2  
Electricity Charging time: 2.5-3 hours  
TPC data cable + manual + 4 massage heads Other



Small round head/for joint

Small fat head/for all parts  
of the body

U head/for vertebrae

Ball head/for large  
muscle groups

### instructions

1. Power on/off button. (long press to start/long press to turn off)
2. Press the key gear to change.
3. Detachable and replaceable vibrator. When disassembling the vibrator; just pull it out with force.

4. Please select 5V charging equipment for charging.

5. For the first time, please charge for 3 hours. Please connect the data cable to the 5V adapter for charging.

6. Button type: when charging, 6 lights will flash, when full, 6 lights will be on, and when working, the lights will indicate the gear position.

7. Touch type: when charging, the nixie tube displays the electric quantity, while working, it displays

# MAINTENANCE

1. Wipe the surface of the massage gun with a slightly damp towel. Dry the massager with soft dry cloth.
2. When stored for extended periods of time, recharge the battery completely. Keep the massager stored in its original packaging.
3. The massager must be returned off for at least a half an hour after an hour's use. Let the massager cool down for at least 30 minutes before using the machine again.



SAVE THESE INSTRUCTIONS



RoHS  
COMPLIANT



Made in China



Sport Scandinavia Aps.  
Niels Bohrs Vej 2. 9900 Frederikshavn. DK  
Tel: +45 70 126 126 . E-mail: info@spsca.dk